

Nona

INDIBA-one

Face&Body Shaping system

Welcome to the new era

Personalized facial and body care

Discover an unparalleled level of personalization

INDIBA-one Face&Body Shaping system uses 448K HZ RF frequency for deep dialing
Energy penetration, taking into account various clinical indications, skin conditions and body parts,
Provide your patients with unique treatments



Applicable

Lifting and firming skin

Eliminating eye wrinkles

Eliminating dark circles and bags under the eyes

Eliminating facial puffiness

Whitening and rejuvenating skin , tightening pores



Applicable

Abdomen

Dissolving the fat, reducing it's size by 2-6 cm, improving abdominal distension, improving cellulite tissue, improving uterine cold and dispelling body cold, improving menstrual irregularities and dysmenorrhea

Shoulder and neck

Unblocking neck meridians, relieving neck stiffness and soreness, promoting lymphatic drainage

Chest

Lifting sagging breasts, clearing underarm lymph, preventing breast hyperplasia and cancer etc .

Buttocks


Lifting the sagging hips

Limbs

Eliminating sagging inside the upper arm and improving skin relaxation

Back and waist

Eliminating muscle soreness, numbness, swelling and pain, dredging back meridian, Eliminating moisture and cold from the body, eliminating soft tissue adhesion, dissolving fat and shaping it .



INDIBA-one integrates CET&RET RF electric heating technology and CVC technology

The perfect combination has brought about a revolutionary change in facial and body treatment.


What is CET&RET RF technology?

The human body is a current conductor. When CET&RET current flows through the human body, it will cause the physiological effects of the skin machine and muscle tissue, which will produce different therapeutic effects according to the depth.


1) When CET&RET RF on the dermis, it mainly acts on collagen fibers and elastic fibers to stimulate collagen to contract, activate collagen cells, repair broken collagen fibers and elastic fibers, and stimulate the secretion of more new collagen and collagen. So as to achieve firming of the skin, improve wrinkles, and enhance the chin contour line to shape facial contours

2) When CET&RET RF on the tendon layer, it mainly acts on the muscle tissue layer. According to the biological effect of the electric heating technology, it acts on the human muscles to increase the temperature of the deep tissues, promote the expansion of blood vessels and accelerate the circulation of lymph, and accelerate the metabolism. The stagnant water and decomposition products in the tissue area are smoothed and metabolized, so as to achieve the effects of relaxing muscles, anti-inflammatory and reducing swelling.

3) When CET&RET RF on the fat layer, the stubborn cellulite is gradually dissolved through the high-frequency diathermy characteristics of the electric heating technology, so as to achieve the effect of eliminating fat, slimming and shaping. Women with dysmenorrhea and congestion have the effect of dispelling blood stasis and dispelling the formation of the body.



Collagen repair schematic



before

after


INDIBA-one integrates CET&RET RF electric heating technology and CVC technology

The perfect combination makes facial anti-aging and body care and shaping a possibility

What is CVC negative pressure technology?

The pressure that the human body normally bears in the air is atmospheric pressure, and the concept of negative air pressure is the state where the air pressure is less than one atmospheric pressure.

1. Specially designed negative air pressure suction heads are used for different depths of kneading massage stimulation for the skin's epidermis, vascular layer, fat layer and nervous system layer.
2. The kneading effect produced by vacuum air pressure helps increase skin tissue elasticity, fight stubborn fat, effectively decompose and break up stubborn cellulite, convert fat into fatty acids and glycerol, which are excreted through lymphatic system metabolism, and at the same time strengthen tissue metabolism.



Your journey with INDIBA-one has just begun.

We hope that the dealer's business will develop with Beautylight. This is why we have assembled a support package and we provide customized services for our dealers so that our dealers can get the most benefits. Start with the first conversation and be with you as your business starts to see benefits. We make sure you have all the training and excellent after-sales service you need. Because we believe that everything is possible!

1

Buy demo test

We will give distributors a very good sample price and buy the first machine to test various performances of the machine.

2

Sign OEM contract

When the customer is satisfied with the test, we will sign a cooperation agreement according to different countries and formally reach an OEM agreement

3

Design new shell

We will design a brand-new machine shell for the dealer to ensure that this shell belongs only to the dealer and keep commercial secrets

6

Form long-term cooperation

After completing the above process, we will establish a long-term cooperative relationship with customers, and we are committed to providing customers with the best after-sales service and technical support

5

Test again

After the customer receives the customized sample, test again to ensure the best satisfaction

4

Send demo to agent

After making the customized sample machine, it will be displayed and tested online. After the customer is approved, it will be sent to the dealer for testing

System Specifications

RF function	CET&RET
Power	300W
RF frequency	448khz
RF Treatment headpiece material	ceramic
Ceramic headpiece size	Ø80、Ø70、Ø50、Ø30
CET%RET headpiece material	medical stainless steel
CVC headpiece size	Ø80、Ø70、Ø60、Ø50
Treatment time	adjustable
Negative pressure	adjustable
Negative pressure frequency	adjustable
Treatment feedback system	open
Rental function	open
Machine size	355*375*1250 (mm)
Work frequency	50-60 Hz
Input voltage	200--240V

Warnings and Risks: IDIBA-one unsuitable for anyone who is pregnant, has an infection or illness, or is heat-sensitive. Cuts, wounds, piercings and tattoos must be avoided and areas with injectable fillers and toxins left for a month before treatment. Patients with a metal implant or implanted device (e.g pacemaker) must consult with a qualified physician first. If deemed suitable, the return pad must be placed Insulate patients from metal objects, especially any in direct contact (e.g. the bed). consultation is always advisable; for contraindications and risks in full, see the IDIBA-one manual